

Use of the ThermoSuit® System for Treatment of Severe Hyperthermia

Robert B. Schock, Ph.D.

From the Department of Research and Development, The Sid Wolvek Research Center,
[Life Recovery Systems](#) HD LLC, Kimmelon, NJ USA

Current medical guidelines recommend cold water immersion as the preferred treatment for heat stroke. The ThermoSuit System is the only FDA-cleared patient cooling device that uses this approach to cool the patient. Independent clinical testing has verified that the ThermoSuit System provides the most rapid cooling of victims of exertional hyperthermia. [LRS Technical Report, April 2017).

Background

Severe hyperthermia is increasingly recognized as a critical health issue. Record high temperatures have been reported in many parts of the world in recent years, and these in turn have been associated with significant increases in rates of death. Thousands of people die each year due to heatstroke, according to the World Health Organization¹. Military personnel operating in high-temperature environments are at a particularly high risk for heat-related disability, injury, or death^{2,3}. Other situations in which high risks of heat injuries are present include the annual Hajj pilgrimage in Saudi Arabia^{4,5} and “rave” dance festivals, in which high heat conditions and recreational drugs such as MDMA have proven to be a lethal combination⁶.

MDMA (otherwise known as “ecstasy”), enhances heat production, induces vasoconstriction, and worsens the body’s ability to survive high environmental temperatures. MDMA can induce fatal hyperthermia in a moderately warm environment. The drug promotes water retention and brain edema, which can be worsened by excessive intake of fluids⁷. The common wisdom that dictates additional fluid consumption during hyperthermia may be ill advised under these conditions.

Heatstroke is a condition resulting from prolonged exposure to intense heat, characterized by high fever and in severe cases convulsions, coma, and death. It is a medical emergency requiring immediate measures to rapidly reduce core temperature. It has been reported that immersion in ice water is superior to all other cooling methods for treatment of severe hyperthermia, including i.v. fluids, cold air, and ice packs^{8,9,10}. Research conducted by a Canadian military research group showed that cold water immersion was particularly effective when the water temperature was

at 2°C, a temperature at which shivering was seldom observed¹¹.

Current Guidelines

Current evidence-based medical guidelines, issued in 2014¹², recommend the use of cold water immersion as the preferred treatment of heat stroke.

The expert panel which developed these guidelines gave cold water immersion a *Recommendation Grade 1A* rating, for treatment of patients (the only cooling treatment to receive the 1A rating). Other cooling methods all received inferior ratings (see Table 1 below)


Treatment Method	Recommendation Grade
Cold-water immersion	1A
Evaporative and Convective Cooling	1C
Cold intravenous fluids,	1C
Body cavity lavage or intravascular cooling devices	2C
Pharmacological Treatment	1B

The cold water immersion method, while highly effective, has until recently not been practical for use while performing other emergency medical treatments such as cardiopulmonary resuscitation. The introduction of the LRS ThermoSuit System has changed that.

Cooling Methods

The ThermoSuit System is the only product that provides controlled patient cooling by ice water immersion – the preferred cooling method recommended in current treatment guidelines, and the only one receiving a 1A recommendation. It delivers a rapidly flowing thin layer of ice water directly over the body in a contained system, rapidly lowering the body's core temperature in minutes. It is designed to be used on a supine patient on a stretcher, facilitating the delivery of CPR and other medical treatments and monitoring procedures simultaneously with cooling. The water is above the freezing point, and thus does not introduce the risk of frostbite that ice-based cooling methods carry.

A comparison of prior studies^{9,13} shows that **ice water immersion can cool victims of severe hyperthermia 40 times faster than standard cooling blankets.** Tan et al¹⁴ recently conducted a clinical study in which hyperthermic subjects were cooled by three cooling methods, including the ThermoSuit System. They demonstrated that the ThermoSuit produced significantly faster core body cooling than did the BCU (Body Cooling Unit) and EMCOOLS devices. The authors cited McDermott et al¹⁰, who stated that ideal cooling modalities provide a cooling rate equal to or greater than $0.155^{\circ}\text{C}/\text{min}$. **Only the ThermoSuit was able to provide a sustained cooling in the ideal range:**



Case Examples

Clinical experience confirms the safe and powerful cooling capabilities of the ThermoSuit system in the treatment of severe, life-threatening hyperthermia. Listed below are several case examples.

Case 1 - Severe Heatstroke: A 47 year old male was found unconscious on the side of a road at midday in Louisiana with the ambient solar exposure + heat index estimated at 101°F. His body temperature was 108 °F. He was brought to a nearby hospital. After attempts at surface cooling produced marginal results, the patient was placed in the ThermoSuit System (TSS) with a starting core temperature of 105° F and cooled to 90.8 °F in 26 minutes, then removed from the TSS and gradually rewarmed. The patient's impending acute renal failure, rhabdomyolysis, and coma completely resolved. Because of elevated cardiac enzymes he underwent cardiac catheterization which was normal and the patient was discharged fully recovered on hospital day five.

Case 2 - Malignant Hyperthermia: A male patient in his 30's developed malignant hyperthermia during surgery (108°F). The anesthesiologist requested that patient be cooled with ThermoSuit. The patient was rapidly cooled to normothermia, and was later discharged with full recovery.

Case 3 - Severe Fever: A female patient in her 50's with non-hodgkins lymphoma arrived in the emergency department with a fever of 105°F, with altered mental status, and hemodynamically unstable. Physicians were unable to control her fever as it rapidly climbed to 107°F. A physician elected to deploy the ThermoSuit. Upon treatment the patient temperature dropped from 106.7°F to 96.8°F in 45 minutes. The next day the patient was awake and was extubated. The patient recovered rapidly, and was later discharged with excellent outcome and full neurologic recovery.

Case 4 - Fever (obese patient): A 44 year-old 298 lb female patient developed a 107.7°F fever. Obese patients tend to be difficult to cool due to the large body mass in relation to surface area. She was cooled with the ThermoSuit to 100° F in only 30 minutes. After being removed from the ThermoSuit, her temperature continued to decrease to 98°F in 5 more minutes.

Discussion

The heatstroke patient described in case 1 above was treated to a level of therapeutic hypothermia (TH) instead of normothermia. His rapid and complete recovery suggests that thermal cell injury may respond to TH in similar manner to ischemic injury and deserves scientific study. It is well known that rapid reversal of core hyperthermia will stop thermal injury to the cells.

The ThermoSuit System uses circulating cold water to cool approximately 90% of the body's surface area, giving it an advantage in terms of heat transfer area vs. all other cooling methods. The ThermoSuit System operates with water that is colder than 10°C (the ThermoSuit usually cools with circulating water at about 2°C). Clinical research by Proulx et al¹¹ demonstrated that colder water significantly increased the cooling rate of immersed volunteers. This is in part due to the increased thermodynamic advantage of using colder water, but a physiological effect is also a factor: Proulx et al reported that only one of seven volunteers shivered when immersed in 2°C water, while six of seven shivered in 8°C water. Thus, the colder water suppressed the shivering response in most subjects. This study reported that subjects in 2°C water lost heat about 50% more rapidly than those in 8°C water.

Shivering raises metabolic rate significantly and adds to the stress on the body. If the goal of cooling is to bring the core temperature to a safe level while minimizing the stress of the cooling process, a rapid cooling induction is desirable. Rapid cooling has also been shown to be desirable when therapeutic hypothermia is indicated^{15,16,17,18,19,20}.



The LRS ThermoSuit® System: The only patient cooling device that provides cold water immersion therapy. This produces the highest noninvasive patient cooling power available.

In summary, the ThermoSuit System has a cooling power greater than any other noninvasive cooling method. This is due to a highly efficient liquid convection heat exchange mechanism that maximizes thermal transfer while minimizing shivering. It uses liquid above the freezing point and does not introduce a risk of frostbite. It enables simultaneous use of medical treatments such as CPR. The ThermoSuit System should be considered as the standard of care for treatment of severely hyperthermic patients.

FDA-Cleared Indications for The LRS ThermoSuit System:

Temperature reduction in patients where clinically indicated, e.g. in hyperthermic patients.

References

1

<http://www.who.int/mediacentre/factsheets/fs266/en/>

2

<http://www.telegraph.co.uk/news/uknews/defence/10178631/Two-soldiers-die-in-heat-while-training-in-Brecon-Beacons.html>

3

<http://phc.amedd.army.mil/topics/discond/hipss/Pages/HeatinjuryPrevention.aspx>

⁴ Al-Ghamdi SM1, Akbar HO, Qari YA, Fathaldin OA, Al-Rashed RS, "Pattern of admission to hospitals during muslim pilgrimage (Hajj)", *Saudi Med J.* 2003 Oct;24(10):1073-6.

⁵ Ganjeh M, Einollahi B, "Mass Fatalities in Hajj in 2015", *Trauma Mon.* 2016 November; 21(5):e43253

⁶ Lin RG, Parker K, "These are the stories of 29 rave-goers who died of drug-related causes", *Los Angeles Times, Dec. 1, 2016.*

⁷ Kiyatkin EA et al, *The Journal of Neuroscience* (2014); 34(23): 7754-7762.

⁸ Hadad E, "Heat Stroke – A Review of Cooling Methods", *Sports Med* 2004; 34(8): pp. 501-511.

⁹ Casa DJ, McDermott EC, Le SW et al, "Cold Water Immersion: The Gold Standard for Exertional Heatstroke Treatment", *Exerc Sport Sci Rev* 2007, 35; 3:141-149.

¹⁰ McDermott BP, Casa DJ, Ganio MS et al, "Acute Whole-Body Cooling for Exercise-Induced Hyperthermia: A Systematic Review", *J Athl Train.* 2009; 44(1): 84-93.

¹¹ C. I. Proulx, M. B. Ducharme, and G. P. Kenny, "Effect of water temperature on cooling efficiency during hyperthermia in humans", *J Appl Physiol* 2003 94(4):1317-1323.

¹² Lipman GS et al, "Wilderness Medical Society Practice Guidelines for the Prevention and Treatment of Heat-related Illness: 2014 Update", *Wilderness & Environmental Medicine* 2014; 25: S55-S65.

¹³ Hart LE et al, "Exertional heat stroke: the runner's nemesis", *Can Med Assoc J* 1980;122(10): 1146-1150.

¹⁴ Tan PMS, Eunice Y. N. Teo EYN, Ali NB, Ang BCH, Iskandar I, Law LYL, Lee JKW, "Evaluation of Various Cooling Systems After Exercise-Induced Hyperthermia", *Journal of Athletic Training: February 2017*, Vol. 52, No. 2, pp. 108-116.

¹⁵ Wolff B et al, "Early achievement of mild therapeutic hypothermia and the neurologic outcome after cardiac arrest", *International Journal of Cardiology* 2009.

¹⁶ Mooney MR, Unger BT, Boland LL et al, "Therapeutic Hypothermia After Out-of-Hospital Cardiac Arrest: Evaluation of a Regional System to Increase Access to Cooling", *Circulation* 2011, 124:206-214.

¹⁷ Sendelbach S et al, "Effects of Variation in Temperature Management on Cerebral Performance Category Scores in Patients Who Received Therapeutic Hypothermia Post Cardiac Arrest". *Resuscitation* 83 (2012), 829-834.

¹⁸ Polderman K and Herold I, "Therapeutic hypothermia and controlled normothermia In the intensive care unit: Practical considerations, side effects, and cooling methods", *Crit Care Med* 2009.

¹⁹ Howes D et al, "Rapid Induction of Therapeutic Hypothermia using Convective-Immersion Surface Cooling: Safety, Efficacy, and Outcomes", *Resuscitation* 2010 Vol. 81, no .4, 388-392.

²⁰ Schock Robert B., Janata Andreas, Peacock W. Frank, Deal Nathan S., Kalra Sarathi, and Sterz Fritz, "Time to Cooling Is Associated with Resuscitation Outcomes", *Therapeutic Hypothermia and Temperature Management.* October 2016, ahead of print. doi:10.1089/ther.2016.0026.

Life Recovery Systems
170 Kinnelon Road
Kinnelon, NJ 07405
Phone: 973-283-2800
Email: info@life-recovery.com
Website: www.life-recovery.com

PN 53164, Rev. C